



## ***EZ-Advancer®***

*Our patented MTG EZ-Advancer® locking valve prevents retraction of the catheter during insertion. This ready-to-use closed system is both easy-to-use and UTI protective.*

*For users on the go, our Mini-Pak version fits in pocket or purse. A privacy bag is included for discreet disposal.*

## MTG EZ-Advancer®

REF #	FR	DESCRIPTION	CASE QTY
31112	12	16" catheter length, firm, with kit components (PVP swabs)	100
31114	14	16" catheter length, firm, with kit components (PVP swabs)	100
31116	16	16" catheter length, firm, with kit components (PVP swabs)	100
31212	12	16" catheter length, soft, with kit components (PVP swabs)	100
31214	14	16" catheter length, soft, with kit components (PVP swabs)	100
32112	12	16" catheter length, firm, with kit components (BZK swabs)	100
32114	14	16" catheter length, firm, with kit components (BZK swabs)	100
32116	16	16" catheter length, firm, with kit components (BZK swabs)	100
32212	12	16" catheter length, soft, with kit components (BZK swabs)	100
32214	14	16" catheter length, soft, with kit components (BZK swabs)	100
30112	12	16" catheter length, firm, without kit	100
30114	14	16" catheter length, firm, without kit	100
30116	16	16" catheter length, firm, without kit	100
30212	12	16" catheter length, soft, without kit	100
30214	14	16" catheter length, soft, without kit	100


## MTG EZ-Advancer® Mini-Pak

REF #	FR	DESCRIPTION	CASE QTY
42112	12	16" catheter length, firm, with kit components (BZK wipe)	100
42114	14	16" catheter length, firm, with kit components (BZK wipe)	100
42116	16	16" catheter length, firm, with kit components (BZK wipe)	100
40112	12	16" catheter length, firm, without kit	100
40114	14	16" catheter length, firm, without kit	100
40116	16	16" catheter length, firm, without kit	100

Full kit contains sterile MTG EZ-Advancer® catheter, pair of latex-free gloves, underpad, gauze pad, BZK or PVP cleansing prep, and privacy bag. Privacy bag is included only in Mini-Pak versions.

## FEATURES:

- The catheter is pre-lubricated, so setup time is minimal.
- The introducer tip is designed to reduce the risk of recurrent Urinary Tract Infections.
- The EZ-Advancer® valve secures the catheter tube after each forward advancement, preventing retraction during insertion.
- Soft, fire polished, smooth eyelets (drainage holes) allow for comfortable passage through the urethra and into the bladder.
- Offset eyes, located near the tip of the catheter provide drainage from two points within the bladder.
- The catheter bag protects the sterility of the tube and collects the urine.
- Access finger holes in the collection bag make holding during filling and draining easier.
- Sterile sample can be collected. The bag capacity is 1500 ml.
- The "Tear here to empty" perforation is easy to rip, allowing for fluid to be poured into the toilet.
- The catheter is available in 12, 14 and 16 Fr., 16" length (see MTG Cath-Lean® should shorter, female length be desired).
- 100% latex-free, DEHP-free catheter tube
- Mini-Pak outer package for travel includes privacy bag for discreet disposal of used supplies



**Medical Technologies of Georgia, Inc.**  
 15151 Prater Drive, Covington, Georgia 30014  
 Tel: 888-511-4239 www.medtechga.com